



FACE SHIELD

PRODUCT CODE: TEC-CAPVISOR _1

DISPOSABLE FACE SHIELD FOR BASEBALL CAPS

FACE SHIELD FOR BASEBALL CAPS

Tecman's lightweight face shield designed for baseball caps provide protection for the face and eyes from splashing or spraying of infectious or hazardous substances. The universal visor slot fits perfectly over the curved visor on baseball caps and is held in place with a soft elasticated strap.

An area is available on the sides of the visor for customised branding and writing names.

The face shields are manufactured from optically clear polyester with an anti-fog coating for maximum visibility. The soft elasticated strap gently curves the clear face shield around the face to provide both maximum shielding and long-term comfort for the wearer. The product can be worn with or without additional eyewear or face protection, such as prescription or protective eyewear and face masks.

Our registered design face shield is manufactured using a patent-protected process.



KEY FEATURES

- Designed to fit perfectly over the visor of baseball caps
- Customised branding and writable nameplate area
- Anti-fog protection
- Optically clear visor for maximum visibility
- Soft elasticated strap for long term comfort
- Lightweight design protects the face and eyes from splashing or spraying
- Can be worn with or without eyewear protection and face masks
- CE Category I



SPECIFICATIONS

SPEC	DETAIL
PRODUCT COLOUR	OPTICALLY CLEAR (NO TINT)
PRODUCT TYPE	FACE SHIELD FOR BASEBALL CAPS
SIZE	330MM X 220MM
OPTICAL CLASS	1 (+/-0.06)
MECHANICAL STRENGTH	S
FIELD OF USE	DROPLETS OR SPLASH OF LIQUIDS
APPLICATION	MEDICAL, GENERAL
STANDARDS	BS EN166:2002 CATEGORY I – SIMPLE DESIGN
CUSTOM BRANDING AREA	100MM X 25MM ON BOTH SIDES OF THE VISOR
SHELF LIFE	2 YEARS
STORAGE CONDITIONS	STORE IN ORIGINAL PACKAGING OUT OF DIRECT SUNLIGHT AT 23°C AND 50% RH

MATERIAL DATA

MATERIAL	DESCRIPTION	SIZE
VISOR FILM	OPTICAL GRADE PET	0.175MM THICK
ELASTIC HEADBAND	MEDICAL GRADE WOVEN ELASTIC	19MM THICK

INSTRUCTIONS FOR USE

With the cap placed on the head, pull the elastic strap of the cap visor around the back of the cap, pulling the visor in front of the face. Position the peak of the cap through the universal visor slot until firmly fixed in place.

ADDITIONAL SPECIFICATIONS

Medical grade materials / Single use, cleaning not recommended / Latex Free / Can be used for prolonged periods as it does not create distortion or fatigue / Universal visor slot is designed to fit perfectly over peaked visors on baseball caps / Visor suitable for use with prescription lenses and surgical masks / Always check before using / Manufactured in the UK

TECHNICAL DRAWINGS

